


PRIGODNA AKCIJA POLOVNE / DEMO OPREME !

The SL-D series employs 12.5x hi-eyepiece eyepieces and Topcon's converging binocular tubes with a parallel magnification body for remarkably sharp and comfortable stereoscopic observation. A wide assortment of optional accessories further expands the clinical efficiencies of the 'D' series.

SL-D7



~~7.250.- €~~
4.500.- € + PDV ✓

Five magnification digital slit lamp

The SL-D7 incorporates a parallel drum magnification changer that allows for quick selection among five convenient magnifications of 6X, 10X, 16X, 25X and 40X. When the optional FD-21 Flash Device is incorporated, the SL-D7 converts to a "true" digital photographic slit lamp providing exceptional images of the cornea, endothelium and fundus to enhance accurate documentation.

DC -3

- » The DC-3 is the perfect combination for still imaging and live image
- » Integrated imaging solution with a 8.3 Megapixel CCD sensor
- » Equipped with a resolution of 1 Megapixel live image
- » Easy installation plug and play
- » Easy capture with joystick button on slit lamp base
- » The DC-3 comes with free IMAGENet i-base™ basic software

~~3.990.- €~~

~~2.500.- €~~
2.500.- € + PDV ✓

